

NOVA-400 DP

CAUTERY MACHINE Unit 400 W DIGITAL

Nova 400 DP Programmable electro surgery generators gives best price performance ratio in similar class of generators offered by other manufacturers and can be widely used in procedures including Urology, GI, Gynecologic, Orthopedic, Thoracic and Plastic Surgery.

Main Features

- Pure Cut Two Blend Modes gives the surgeon varying degree of Homeostasis in the cut mode. Endo cut mode is provided for added flexibility.
- True Spray Coag mode offers massive and rapid non-contact coagulation. Fulgurate Coag mode for sealing bleeders. Desiccate Coag mode for contact coagulation.
- Micro / Macro Bipolar Modes provides controlled desiccation without charring, adhesion and blanching of tissues with bipolar cut and coag facility. Bipolar cut mode allows bipolar cutting operations.
- Separate Monopolar and Bipolar R.F. output system allow two surgeons to work separately.
- Absolutely Safe as Isolated RF Output provides direct current return to the generator significantly reducing the risk of alternate site burns.

