

SI-117A

ELECTRIC EXAMINATION TABLE

Backed by experience team, we are presenting a wide series of 3 Section Examination Couch. These products are extremely admired by patrons due to their strong nature and low maintenance. Our product is built up employing the best quality material which is procured from trustworthy merchants of market.

Specifications:

- Height adjustment by electromotive drive
- The electric motor is controlled by 1 foot switch (possible option: Hand switch)
- Head part step less up- (50°) and downwards (-45°)adjustable by gas spring.
- Length head part: 450 mm
- Middle part fixed. Length: 1.290 mm
- Foot part step less adjustable from downwards vertical position (-90°) to upwards position (28°) by gas spring.
- Length foot part: 450 mm
- Upholstery base: Laminated wood
- Length min.: 1.790 mm, Length max. 2.190 mm

