



## ICU BED MECHANICAL ABS PANEL COLLAPSIBLE RAILING & MATTRESS

**Code : ASI-107**

### **Description :**

- Frame work made of Rectangular M.S. Tube
- 4 Section Top made of perforated M.S. Sheet
- Back Rest, Knee Rest, TB/RTB Positions maneuvered by Separate screw from Foot End.
- ABS head and foot Panels.
- Collapsible Aluminum safety side Railings.
- SS IV Rod with four Locations
- Bed Mounted on 150mm Dia Castor (2 with Brakes)
- Overall size: 210L x 90W x 60H cms
- Finish: Epoxy Powder Coated
- Optional : Four Sectional Mattress