



LABOUR TABLE IMPROVED MODEL

Code : ASI-136

Description :

- Approx. Dimensions: 72"Lx27"Wx30"H.
- Tubular frame work mounted on PVC stumps.
- Stainless Steel top in three sections. U- cut at foot end.
- Trendlenburg positions adjustable on screw mechanism by handle.
- Backrest section adjustable on ratchet system.
- A pair of lithotomy Rods with straps with height adjustable.
- Complete with self locking traction handle at side & I. V. Rod.
- Fine finished.