



Sodium Combination Electrode

Application :

Steam condensate, meats, fish, diary, fruit juices, fermentation, body fluids, ground/sea water, soils,

Part No	ISE-23
Ion Type	Sodium
Symbol	Na ⁺
Make	LABMAN
Refillable or Gel (Sealed)	Refillable
Concentration range	10 ⁻¹ M to 10 ⁻⁶ M
Limits (ppm)	2300 ppm to 0.02 ppm
pH Range	above 10
Temperature Range	5 to 60 °C
Electrode Resistance	less than 1 megohms
Reproducibility	±2%
Membrane	Solid State
Material	Glass Electrode
Slope	+54 to +60 mV
Cable Length	1 meter
Connector	BNC
Main Interferences	Li ⁺ , K ⁺
Dimension	Φ 12 x 165 mm
Accessories	Sodium Standard Solution & ISAB (500ml each)