



### Salient Features

- Based on E.M.F.C. Technology
- A built in under hook for below weighing
- Density measurement (Optional)
- RS232C Interface
- Tare range upto max. Capacity of the Balance
- High Resolution, Quick Weighing, Accurate Result
- Automatic External Calibration
- Multifunction Weighing Units
- S.S. Weighing Pan, Level Indicator, Adjustable Fit
- Counting Weighing Mode

## High Precision Balance

MODEL	HPB 2000	HPB 3000	HPB 5000
Capacity	2000 g	3000 g	5000 g
Make	WENSAR		
Readability	0.01 g ( 10 mg)		
Repeatability ( ± )	0.01 g		
Linearity ( ± )	± 0.02 g		
Display	VFD ( Vacuum Fluorescent Display )		
Pan Size	150 mm Ø		
Response Time	3 - 4 Sec		
Calibration	Automatic External		
Tare Range	Full		
Operating Temperature	10°C to 40°C		
Interface	RS232C		
Power Supply	AC 220V / 50-60 Hz		
Dimension (LxWxH) mm	530 x 435 x 345		
Weight	6.9 kg	★ p o n a : n e a e , r o e y s a n	

SPECIFICATIONS

★ Specifications are subject to change without notice.